# M.I.C.EVALUATOR™

X7216B

UK

**M.I.C.Evaluator™****Intended Use**

M.I.C.Evaluator™ (M.I.C.E.) is a system for quantitatively determining the Minimum Inhibitory Concentration (MIC) of an antibiotic against a test organism. The Oxoid M.I.C.E. provides a gradient of antibiotic stabilized on a plastic strip with 30 graduations, to give an accurate MIC over the range 256 µg/ml – 0.015 µg/ml (specialist high and low level strips are also available).

**Principle of the Test**

On application of the M.I.C.E. to pre-inoculated agar, the antibiotic starts to release from the plastic forming a defined concentration gradient in the area around the strip. After appropriate incubation of the test, a zone of inhibition will have formed with the M.I.C.E. at the centre. The MIC can be easily read, using the graduated scale, at the point where the growth of the test organism touches the M.I.C.E..

**Components**

M.I.C.E. are supplied in individual peelable sachets with a desiccant to maintain the stability of the antibiotic on the strip. A variety of pack sizes is available including 10 and 50 strips of a single antibiotic or recognised antibiotic combination, in a rigid outer carton with an instruction leaflet.

**Materials required but not provided**

Agar plates with appropriate media, Inoculum suspension medium, Sterile loops and swabs, Sterile forceps, McFarland turbidity standards, Incubator, Modified atmosphere environments for CO<sub>2</sub>, microaerophilic and anaerobic conditions, Quality control strains.

**Precautions**

This product is for *in vitro* diagnostic use only. Do not use the product beyond the stated expiry date. Do not use the strips if they show signs of damage – this may prevent sufficient contact of the strip with the agar and will hinder the formation of the gradient. Do not use if the seal on the sachet is broken. The product contains antibiotic – only handle the strip at the top above the white chevron where the logo and antibiotic code are located (diagram 1a). Standard guidelines for the safe handling and disposal of infectious organisms should be observed throughout all procedures.

**Storage and Handling**

Ensure that the sachets are stored so that they are not damaged or bent. The unopened sachets should be stored at 2 - 8°C. If a sachet is opened the strip should be used within 15 minutes – unused strips should be discarded after this time. **It is recommended that the sachets be allowed to come to room temperature before opening so that condensation does not collect on the strip which may result in a disturbance of the gradient.** The sachets are opened by peeling gently at the corners indicated with the blue arrow, this will reveal the M.I.C.E. inside the sachet. Handle the M.I.C.E. only at the top where the logo and antibiotic code are as shown in diagram 1a – DO NOT TOUCH THE ANTIBIOTIC GRADIENT AREA OF THE STRIP.

**Quality Control Procedures**

It is recommended that control strains are tested at appropriate intervals; this can either be when each test is performed, or as recommended in the guidance from the national reference groups for antibiotic susceptibility testing. The control ranges for various antibiotic/organism combinations are listed in Table 2 or refer to local guidelines. The quality control ranges are based on ISO 20776-1 (2006). If the MIC obtained for a control strain organism against a specific antibiotic is outside the range specified, the strips should not be used until the reason for the error has been determined (media, depth, inoculum level, incubation conditions or deterioration of the strip).

**Method****Media**

Select the appropriate medium for the organism to be tested from Table 1. All media should be poured to a depth of 4mm +/- 0.5mm, and conform to appropriate quality controls before use.

**Inoculum**

Prepare the appropriate McFarland inoculum level indicated in Table 1 by emulsifying several colonies from a pure culture into a suitable suspension medium and compare the turbidity to the appropriate McFarland standard. The suspension should be used within 15 minutes of preparation. Dip a sterile cotton swab into the suspension and remove the excess moisture by pressing against the edge of the tube. Inoculate the plate by swabbing in at least three different directions, taking care to do this evenly to ensure no gaps are left in the deposited inoculum. Allow the surface of the agar to dry completely before applying the M.I.C.E., however the strip should be applied within 15 minutes of inoculation to avoid pre-growth of the organism.

**Application of the M.I.C.E.**

Using sterile forceps, remove the M.I.C.E. from the sachet by handling the end with the logo and antibiotic code. Place the strip with the scale facing upwards, and the antibiotic gradient downwards in contact with the agar. Apply the strip by putting the end with the lowest concentration onto the plate first and then carefully rolling the strip onto the agar to ensure good contact with the entire length of the M.I.C.E.. If air bubbles do become trapped under the strip, gently ease them to the edge using sterile forceps taking care not to move the strip. If the M.I.C.E. is placed in an incorrect position, DO NOT MOVE THE STRIP since antibiotic will be released into the agar immediately. It is recommended that no more than two M.I.C.E. should be applied to a 90mm plate, and no more than six to a 140mm plate. Once the strips are applied to the agar the plates should be incubated under the appropriate conditions immediately to avoid pre-diffusion of the antibiotic.

**Incubation**

Plates should be incubated in an inverted position under the conditions appropriate for the organism as indicated in Table 1.

**Interpretation of Results**

In order to assist with clinical interpretations only the full dilutions are indicated on the M.I.C.E. strips (white sections). The half step dilutions are the black sections, and if required can be determined as the midpoint between the values indicated in the white sections, or refer to [www.oxoid.com](http://www.oxoid.com) for a detailed diagram. If the growth touches the strip on a white section, the MIC is read as the value in that section as shown in diagram 1b. If the growth touches the strip in a black section then for clinical interpretation the MIC is read as the value in the next white section above where the growth has finished (diagram 1c). If there is growth along the entire length of the strip (no zone of inhibition), the MIC should be read greater than the highest value on the strip. If a large zone of inhibition is obtained, and the growth of the organism does not intersect with the strip, the MIC should be read as less than the lowest value on the M.I.C.E.. Detailed information regarding the interpretation of the growth patterns can be obtained at [www.oxoid.com](http://www.oxoid.com).

**Performance Characteristics**

Performance characteristics generated from internal trial data are shown in Table 3

**Limitations**

The M.I.C.E. results are *in vitro* results, and therefore only give an indication of the *in vivo* susceptibility of the test organism. The decision to use an antibiotic for therapy against the test organism is the responsibility of the clinician who will take into account other factors which may influence the *in vivo* activity of the compound.

**Packaging Symbols**

☒	For <i>in vitro</i> diagnostic use	☐	Batch number of the M.I.C.E.
Ⓜ	Expiry date	Ⓜ	Storage temperature
☒	Number of tests per package	☒	Instruction leaflet inside the pack
		☐	Conformité Européene

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ATCC® catalogue marks are trademarks of ATCC \* ranges derived from the CLSI standard
Oxoid Ltd. Wade Road, Basingstoke, Hants, RG24 8PW, England