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1. Croxatto A, Prod'hom G, Faverjon F *et al.* Laboratory automation in clinical bacteriology: what system to choose? *Clin Microbiol Infect.* 2016; 22:217-35.
2. 李珮寧，邱晟榮，林函頤，蔡文城。從 PreLUD 自動化檢體接種機與手工接種痰檢體的初步效能比較探討市面上的三種自動化接種系統。檢驗及品保雜誌。2017; 6:63-70。
3. 邱彥昕，謝賢修，羅晟展等人。評估 CMP™ 嗜氧檢體輸送管對臨床細菌之保存效果。檢驗及品保雜誌。2012; 1:1-6。
4. Dube FS, Kaba M, Whittaker E *et al.* Detection of *Streptococcus pneumoniae* from different types of nasopharyngeal swabs in children. *Plos ONE* 2013; 8:e68097.
5. Dalmaso G, Ferrari M, Paris A. Validation of the new Irradiated Nylon™ flocked QUANTISWAB™ for the quantitative recovery of micro-organisms in critical clean room environments. PDA Europe Conference and Exhibition "Best Practices in Aseptic Manufacturing". 2007. Milan, Italy.
6. Copan. FLOQSwab®: Frequently Asked Questions. 2022. <https://www.copanusa.com/floqswab-frequently-asked-questions/#Q1>
7. 葉芝榕，歐柏廷，蔡文城。尼龍植絨、人造絲與棉花拭子採集微生物檢體的效能。檢驗及品保雜誌。2022; 10:162-6。
8. Clinical and Laboratory Standards Institute (CLSI). Quality control of microbiological transport systems, 2nd ed. Approved standard M40-A2. 2015. CLSI. Wayne, PA, USA.
9. Thermo Fisher Scientific Manual. Amies transport medium w/ and w/o charcoal 2008.
10. Olsen CC, Schwebke JR, Benjamin Jr WH *et al.* Comparison of direct inoculation and Copan transport systems for isolation of *Neisseria gonorrhoeae* from endocervical specimens. *J Clin Microbiol* 1999; 37:3583-5.
11. Lindy G, Frans R, Duduzile N *et al.* Evaluation of the Copan eSwab®, a liquid-based microbiology transport system, for the preservation of *Neisseria gonorrhoeae* at different temperatures. *South Afr J Infect Dis* 2017; 1:1-4.
12. Serra-Pladevall J, Gulin BC, Vila ON *et al.* Preservation of *Neisseria gonorrhoeae*: should swabs be refrigerated or not? *Neisseria gonorrhoeae* preservation. *J Microbiol Methods* 2018; 145:37-9.
13. Wade JJ, Graver MA. Refrigeration does not compromise recovery of *Neisseria gonorrhoeae* from charcoal transport swabs. *Sex Transm Infect* 2005; 81:93-4.

CMP efSwab[®]-Amies Broth is Effective in Preserving both Aerobic and Anaerobic Pathogenic Bacteria

Yi-Hsuan Huang¹, Yu-Jung Chang, Wen-cherng Tsai¹⁻³

¹Super laboratory Ltd., New Taipei City ; ²Institute of Microbiology and Immunology, and ³Institute of food Safety and Health Risk Assessment, National Yang Ming Chiao Tung University, Taipei, Taiwan

Abstract

Automated specimen inoculation devices, such as Inoqula (BD-Kiestra, USA) and WASP (Walk-Away Specimen Processor, Copan, Italy), have become widely used in large microbiology laboratories. As these instruments require liquid medium for inoculation, it is common to use flocced swabs coupled with Amies transport medium for specimen collection. To meet the requirement of automated inoculation, the Taiwan Creative Life Science Co., Ltd. has developed a kit, referred to as CMP efSwab[®]-Amies broth, for specimen collection. In this study, we evaluated its efficacy according to the CLSI Microbiological Transport System Specification M40-A2 (2014) on both aerobic bacteria including *Neisseria gonorrhoeae* (ATCC 43069) and *Haemophilus influenzae* (ATCC 10211) and anaerobic bacteria including *Peptostreptococcus anaerobius* (ATCC 27337) and *Bacteroides fragilis* (ATCC 25285). To simulate the process of specimen transport, each bacterial suspension was adjusted to McFarland 0.5 with a bacterial count of approximately 1.5×10^8 CFU/mL, and 100 μ L each of such suspensions was placed on a flocced swab (efSwab[®]). The inoculated swab was then

inserted into the tube containing Amies broth and stored at 4°C or room temperature (25°C) for 6, 24, and 48 hours. According to the CLSI standard, the decrease in the log value of bacterial count of a specimen should not be greater than 3 and the increase value should not be greater than 1 after being stored at 4°C or room temperature. Results showed that *H. influenzae*, *P. anaerobius*, and *B. fragilis* could be recovered after 48 hours and *N. gonorrhoeae* could be recovered after 24 hours of storage at 4°C or room temperature with a change in the log value of bacterial count less than 3. As these results met the CLSI M40-A2 standard, we conclude that the CMP efSwab[®]-Amies broth specimen collection kit is compatible with automated specimen inoculation instruments.

Keywords: Automated specimen inoculation devices, efSwab[®]-Amies broth, *Neisseria gonorrhoeae*, *Haemophilus influenzae*, *Peptostreptococcus anaerobius*, *Bacteroides fragilis*