

# 大腸桿菌及大腸桿菌群 酵素呈色檢測法


## Chromogenic ECC agar

### Chromogenic ECC agar (NIEA 配方, E237.52B)

- SteriPack pre-poured, 200 ml 瓶裝 (適用傾注平板法)
- Ready to use 55 mm 平板 (搭配濾膜法)

### Chromogenic ECC agar, modified (NIEA 改良版)

- SteriPack pre-poured, 150 ml 瓶裝 (適用傾注平板法)
- Ready to use 90 mm 平板 (適用塗抹平板法)
- Ready to use 55 mm 平板 (搭配濾膜法)
- Powder form (粉末自行滅菌配製)



### Chromogenic ECC agar, modified NIEA 改良版：

- 適用於水質與食物檢體；濾膜法或傾注法
- 菌落形狀完整，顆顆分明，更不須看產氣
- 呈色不模糊，清楚分別 coliform 與 E. coli

### CMP® SteriPack Pre-poured bottled Medium 無菌瓶裝回融培養基：

- GMP 認證工廠製造
- 適用於採取混稀傾注法 (pour plate) 時使用
- 預先凝固，便於保存與即時使用
- 僅需加熱回融，不用秤重分裝滅菌
- 200 ml or 150 ml 包裝，便於傾倒



致力於微生物領域

