



## Applications:

- QC of microbial identification and detection systems and commercial diagnostic assays including biochemical, PCR and molecular technologies
- QC of CLSI® Antimicrobial Susceptibility Testing (AST) and many other common tests
- Verification and Validation
- Proficiency

## Product Details:

- **Live Culture: KWIK-STIK™ 2-Pack**
  - Two KWIK-STIKs of each microorganism
  - Refrigerated storage
  - For *In Vitro* Diagnostic (IVD) use
- **Inactivated Pellets and Swabs**
  - Multiplexed pellets
  - Swabs mimicking patient sample processing
  - Room temperature storage
  - For *In Vitro* Diagnostic (IVD) use

## Highlights:

- **Diagnostic System and Syndrome-Specific QC Sets and Panels** contain the recommended quality control strains for a particular instrument or assay in an easy-to-use format
- **Regulation-Specific QC Sets and Panels** contain the regulation-specific (e.g., CLSI or EUCAST) recommended quality control strains for a particular assay in an easy-to-use format
- One catalog number for fast, easy ordering
- Protects your reputation by providing documentation and assurance that your laboratory can consistently produce accurate results
- Ensures testing procedures and materials are working properly
- Convenient test-ready formats save you time and money
- Technical Support experts available for guidance

## BD MAX™

- 8193 BD MAX™ CT/GC/TV 20-Day Control Panel
- 8202 BD MAX™ Enteric Parasite 20-Day QC Panel
- 8204 BD MAX™ Enteric Parasite Control Panel
- *C. difficile* Control Panel (Inactivated Swab) - *Coming Soon*
- Comprehensive Enteric Bacterial Control Panel - *Coming Soon*
- CT/GC/TV Control Panel (Inactivated Pellet) - *Coming Soon*
- 8179 Enteric Bacterial Organism Set (Live Culture)
- 8191 Extended Enteric Bacterial Verification Panel (Inactivated Pellet)
- 8211 Enteric Viral Control Panel
- 8210 Enteric Viral Verification Panel
- Group B Streptococcus (GBS) Control Panel (Inactivated Swab) - *Coming Soon*
- Methicillin-Resistant *Staphylococcus aureus* (MRSA) Inactivated Swab - *Coming Soon*
- Methicillin-Susceptible *Staphylococcus aureus* (MSSA) Inactivated Swab - *Coming Soon*
- Methicillin-Susceptible *Staphylococcus epidermidis* (MSSE) Inactivated Swab - *Coming Soon*
- 8174 MRSA Organism Set (Live Culture)
- 8178 SA Organism Set (Live Culture)
- 8173 Strep B Organism Set
- 8209 Vaginal Control Panel
- 8208 Vaginal Verification Panel

## Cepheid GeneXpert®

- 8187 Carbapenem-resistant Enterobacteriaceae (CRE) Control Panel (Inactivated Swab)
- 8200 Cepheid Xpert® *C. difficile* Control Panel
- 8188 Cepheid Xpert® CT/NG Control Panel
- 8194 Cepheid Xpert® GBS LB Control Panel
- 8195 Cepheid Xpert® MRSA/MRSA NxG Control Panel
- 8196 Cepheid Xpert® SA Nasal Complete Control Panel
- 8199 Cepheid Xpert® Respiratory Control Panel
- 8171 *C. difficile* QC Set (Live Culture)
- 8190 Enterovirus (EV) Control Panel (Inactivated Swab)
- 8164 GBS QC Set (Live Culture)

## BD Viper™

- 8214 Human papillomavirus (HPV) Control Panel (Inactivated Pellet)
- 8216 Human papillomavirus (HPV) Verification Panel (Inactivated Pellet)

## BioFire® FilmArray®

- 8215 Blood Culture Identification (BCID) Control Panel
- 8201 Blood Culture Identification (BCID) Verification Panel
- 5229P Blood Culture Identification (BCID) Verification Panel (Live Culture)
- 8217 Respiratory (21 Targets) Control Panel

## Luminex VERIGENE®

- 8184 Enteric Pathogens Control Panel (Inactivated Pellet)
- 5226P Gram-Negative Blood Culture Control Panel (Live Culture)
- 8181 Gram-Negative Blood Culture Control Panel (Inactivated Pellet)
- 5225P Gram-Positive Blood Culture Control Panel (Live Culture)
- 8180 Gram-Positive Blood Culture Control Panel (Inactivated Pellet)
- 8198 Respiratory Pathogens (16 Strains) Control Panel (Inactivated Pellet)

- 8219 Group A Streptococcus (GAS) Control Panel
- 8205 Norovirus Control Panel (Inactivated Swab)
- 0360MSSA *Staphylococcus aureus* subsp. *aureus* derived from ATCC® 25923™\* (MSSA Live Culture Positive Control)
- 0158MRSA *Staphylococcus aureus* subsp. *aureus* derived from ATCC® 700699™\* (MRSA Live Culture Positive Control)
- 0371 MSSE *Staphylococcus epidermidis* derived from ATCC® 12228™\* (MSSE Live Culture Negative Control)
- 8189 *Trichomonas vaginalis* (TV) Control Panel (Inactivated Swab)
- 8203 Vancomycin Resistant Enterococcus (VRE) Control Panel (Inactivated Swab)



\*Look for the ATCC Licensed Derivative® Emblem for all products derived from ATCC® cultures. The ATCC Licensed Derivative Emblem, the ATCC Licenses Derivative Work Mark, and the ATCC Catalog Marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.