

高專一性的凝膠法試劑

PYROSTAR™ ES-F 凝膠法
PYROSTAR™ ES-F Gel Clot products



- Multi-test Gel Clot Lysate, 0.015, 0.03, 0.06, 0.125 and 0.25 EU/ml
- Single-test Gel Clot Lysate, 0.015 EU/ml

PYROSTAR™ ES-F Series products are now available through Lonza. For more details and more ordering information please visit Lonza website: www.lonza.com/pyrostarsf



PYROSTAR™ ES-F Series

內毒素凝膠法試劑與標準品

規格與特色

- 內毒素專一性的裂解液，免除(1→3) β-D-glucan造成偽陽性干擾
- 為客戶成本量身打造不同試驗數量及單支凝膠試劑
- 具備五種不同靈敏度，範圍從0.015-0.25 EU/mL
- 適用各種樣本類型，PPC回收率合乎藥典標準

PYROSTAR ES-F 凝膠法裂解液組

PYROSTAR ES-F 凝膠法裂解液組 (5.2 mL)：4瓶 5.2 mL裂解液+1瓶CSE (500 ng/vial)

Multi-Test Lysate with matched CSE	
Catalog No.	Gel-Clot 靈敏度(EU/mL)
WPEK4-50015	0.015
WPEK4-50003	0.03
WPEK4-50006	0.06
WPEK4-50125	0.125
WPEK4-50025	0.25

PYROSTAR ES-F 單支凝膠試劑

PYROSTAR ES-F 單支凝膠試劑：25支單支凝膠試劑 +1瓶CSE (500 ng/vial)

Single-Test Lysate with matched CSE	
Catalog No.	Gel-Clot 靈敏度(EU/mL)
WPESK-0015	0.015

PYROSTAR 凝膠法內毒素標準品組

PYROSTAR 凝膠法內毒素標準品組(CSE)

Control Standard Endotoxin	
Catalog No.	規格
CSE4037-5006	6瓶 (500 ng/vial)

