


CMP[®] SteriPack Media

Ideal for the environmental monitoring in the cleanroom.

Features:

- Gamma irradiation sterilization (Integrated irradiation color indicator)
- Growth promotion test meets the requirement of USP/EP
- Easy to open
- Add neutralizer to avoid false negative results
- Suitable for pharmaceutical industry

Steps of unpacking:



Item Description	Product description	Specification
SteriPack TSA, Triple	CMP0800011	90 mm Mono, 40 PC/BOX
	CMP0800012	90 mm Mono, 100 PC/BOX
SteriPack TSA w/Lec. and P80, Triple	CMP0800192	90 mm Mono, 100 PC/BOX
SteriPack TSA w/Lec., P80, and Lac. 1U, Triple	CMP0800212	90 mm Mono, 100 PC/BOX
SteriPack TSA w/Lec., P80, and Penase, Triple	CMP0800182	90 mm Mono, 100 PC/BOX
SteriPack TSA with 4 Neutralizers, Triple(Lec., P80, His, Na2S2O3)	CMP0800142	90 mm Mono, 100 PC/BOX
SteriPack SDA, Triple	CMP0800051	90 mm Mono, 40 PC/BOX
Contact Media		
SteriPack Contact TSA w/Lec. and P80, Triple	CMP0800041	55 mm Contact, 40 PC/BOX
	CMP0800042	55 mm Contact, 100 PC/BOX
SteriPack Contact TSA w/Lec., P80, and Lac. 1U, Triple	CMP0800092	55 mm Contact, 100 PC/BOX
SteriPack Contact TSA w/Lec., P80, and Penase, Triple	CMP0800072	55 mm Contact, 100 PC/BOX
SteriPack Contact TSA with 4 Neutralizers, Triple (Lec., P80, His, Na2S2O3)	CMP0800112	55 mm Contact, 100 PC/BOX
SteriPack Contact TSA, Triple	CMP0800021	55 mm Contact, 40 PC/BOX

*Lec: Lecithin, P80: Polysorbate 80, Lac: Lactamator

C20221110